



# RESPIRATORY & ANESTHESIA

**GC MEDICA ENTERPRISE LTD. ( WUXI )**  
GC MEDICA ΕΠΙΧΕΙΡΗΣΗ ΓΙΔ ( ΜΟΥΧΙ )

# CONTENTS

## 01

### OXYGEN THERAPY

● Capnography Mask	01
● CO <sub>2</sub> Nasal Oxygen Cannula	02
● Elastic Head Strap Cannula	03
● Universal Green Oxygen Bubble Tubing	04
● Oxygen Mask	05
● Non-Rebreather Mask	06
● Multi-venturi Mask	07
● Adjustable Venturi Mask	08
● Aerosol Masks with Nebulizer	09-10
● Mouthpiece with Nebulizer	10
● Tracheostomy Venturi Mask	11
● Tracheostomy Mask	12

## 02

### AIRWAY MANAGEMENT

● Berman Airway	13
● Guedel Airway	14
● Nasopharyngeal Airway	15
● Robertazzi Nasopharyngeal Airway	15
● Endotracheal Tube	16
● Reinforced Endotracheal Tube	17
● Suction Plus Endotracheal Tube	18
● Oral Preformed Endotracheal Tube	19
● Nasal Preformed Endotracheal Tube	20
● Right/Left-sided Endotracheal Tube	21
● CELESTIAL Adjustable Tracheostomy Tube	22
● CELESTIAL Tracheostomy tube with NeoSuction	23
● CELESTIAL Reinforced Tracheostomy tube	24
● SERENITY Tracheostomy Tube	25
● SERENITY Tracheostomy Tube Kit	26
● SERENITY Tracheostomy Tube with NeoSuction	27
● SERENITY Tracheostomy Tube Kit with NeoSuction	27
● GRASMERE Fenestrated Tracheostomy Tube	28
● GRASMERE Fenestrated Tracheostomy Tube Kit	28
● AZURE Non-adjustable Tracheostomy Tube	29
● Intubating Stylet	29
● Endotracheal Tube Introducer	30
● Disposable Laryngeal Mask Airway	31
● 90° Preformed Laryngeal Mask Airway	32
● Reinforced Laryngeal Mask Airway	33
● Reusable Laryngeal Mask Airway	34
● Double-lumen Silicone Laryngeal Mask Airway	35

# 03

## RESPIRATORY CARE

- Mucus Specimen Trap with Vacuum Breaker ----- 36
- Mucus Trap with Mouthpiece ----- 37
- Mucus Trap with Tapered Connector ----- 38
- Mucus Extractor ----- 39-40
- Mucus Specimen Trap ----- 41
- Mucus Specimen Traps ----- 42
- Suction Catheter ----- 43-44
- MDI Spacer ----- 45

# 04

## ANESTHESIA

- Corrugated Breathing Circuit ----- 46
- Expandable Breathing Circuit ----- 47
- Coaxial Breathing Circuit ----- 48
- Smoothbore Breathing Circuit ----- 49
- Cather Mount ----- 50
- Anesthesia Mask ----- 51
- Anes-Air Cushion Mask ----- 51
- Eco-Anes Mask ----- 51
- Manual Resuscitator ----- 52
- Incentive Breathing Exerciser ----- 53
- Volumetric Breathing Exerciser ----- 54
- Gas Sampling Line ----- 54
- Heat Moisture Exchange (HME) ----- 55
- Heat Moisture Exchange Filter (HMEF) ----- 56
- Bacterial/Viral Filter ----- 57

# 05

## ABC MOUTHPIECE AND FILTER KIT ----- 58



# 01 OXYGEN THERAPY

## Capnography Mask

- Capnography Mask is a soft pliable contoured design that provides a comfortable fit and incorporates a kink resistant oxygen line, plus CO<sub>2</sub> sample line outlet
- A CO<sub>2</sub> in-line filter option is also available



Product code	Gas sampling line	Size
GC8060A01	without	Adult elongated
GC8060A02	without	Adult standard
GC8060A03	without	Pediatric elongated
GC8060A04	without	Pediatric standard
GC8060B01	with	Adult elongated
GC8060B02	with	Adult standard
GC8060B03	with	Pediatric elongated
GC8060B04	with	Pediatric standard

# CO<sub>2</sub> Nasal Oxygen Cannula

· Adult low-flow plus CO<sub>2</sub> sampling cannula with 2m gas sampling tube for CO<sub>2</sub>



Product code

GC8022L01


Size

Adult

## Elastic Head Strap Cannula

Elastic Head Strap Cannula is designed for a more comfortable feel for those who need daily respiratory therapy. These soft, Over-the-Ear Cannulas are simply inserted into the nostrils to allow proper air delivery from the an oxygen concentrator machine.

- Lightweight, soft, transparent vinyl
- Adjustable elastic band for comfortable, snug fit below the ears
- Available with 7 ft, oxygen supply tubing



Product code

GC8022S11

Size

Adult

# Universal Green Oxygen Bubble Tubing


Universal Green Oxygen Bubble Tubing is a 30 metre length of green tubing designed for clinicians who can cut the tubing to the required length and diameter so that the tubing is matched to the appropriate mask/nebulizer and corresponding oxygen source.



Product code	Description	Length
GC1301A03	OD4.6mm, ID3mm, green	30m

# Oxygen Mask


- A full range is available to meet all requirements in pediatric standard, pediatric elongated, adult standard and adult elongated
- Available with swivel oxygen adapter for easy application and patient's comfort
- Smooth featured edges are designed for anatomical compatibility



Product code	Size	Tubing length	Nose clip	Swivel oxygen connector
GC8015A01	Adult elongated	2.1m	✓	
GC8015A02	Adult standard	2.1m	✓	
GC8015A03	Pediatric elongated	2.1m	✓	
GC8015A04	Pediatric standard	2.1m	✓	
GC8026C01	Adult	2.1m		✓
GC8026C02	Pediatric	2.1m		✓

# Non-Rebreather Mask

- Low resistance, Latex-free check valve allows exhaled gas to escape, helping patient receive better quality air
- A full range is available to meet all requirements in pediatric standard, pediatric elongated, adult standard and adult elongated




Product code	Size	Tubing length	Nose clip	Bag volume
GC802501	Adult elongated	2.1m	✓	1.0L
GC802502	Adult standard	2.1m	✓	1.0L
GC802503	Pediatric elongated	2.1m	✓	0.75L
GC802504	Pediatric standard	2.1m	✓	0.75L



## Multi-Venturi Mask

- Variable sized masks
- 6 color-coded diluters for different concentration
- Includes adapter for high-humidity entrainment

Blue: 24%      Green: 35%  
 Yellow: 28%    Pink: 40%  
 White: 31%     Orange: 50%




Product code	Size	Tubing length	Nose clip
GC8026B01	Adult elongated	2.1m	✓
GC8026B02	Adult standard	2.1m	✓
GC8026B03	Pediatric elongated	2.1m	✓
GC8026B04	Pediatric standard	2.1m	✓



## Adjustable Venturi Mask


- Safe and simple delivery for variable concentration of oxygen
- Color coded, low and medium concentration diluters, green diluters for 24%, 26%, 28%, 30%, white diluter for 35%, 40%, 50%
- Locking ring secures flow settings
- Includes adapter for high-humidity entrainment



Product code	Size	Tubing length	Nose clip
GC8026A01	Adult elongated	2.1m	√
GC8026A02	Adult standard	2.1m	√
GC8026A03	Pediatric elongated	2.1m	√
GC8026A04	Pediatric standard	2.1m	√


# Aerosol Masks with Nebulizer

- Nebulizer is 6cc, designed to be used with any standard nebulizer kit
- Standard fitting allows quick connection to nebulizer kit. Anatomically formed for maximum comfort
- Direct flow design maximizes flow rate and delivers more medication to patient
- Clear transparent shell for easy observation and monitoring



GC802731

Product code	Size	Nebulizer	Tubing length	Nose clip
GC802731	Adult elongated	✓	2.1m	✓
GC802732	Adult standard	✓	2.1m	✓
GC802733	Pediatric elongated	✓	2.1m	✓
GC802734	Pediatric standard	✓	2.1m	✓



GC803401

Product code	Size	Nebulizer	Tubing length	Nose clip
GC803401	Adult elongated			✓
GC803402	Adult			✓
GC803403	Pediatric elongated			✓
GC803404	Pediatric			✓

## Mouthpiece with Nebulizer

· Available with 6cc and 20cc nebulizer




Product code	Nebulizer	Tubing length	Aerosol tube
GC802801	6cc	2.1m	✓
GC802801	20cc	2.1m	✓



## Tracheostomy Venturi Mask







- Specialized venturi style mask allows for the delivery of oxygen to patient who has a tracheostomy
- Equipped with 6" flex tube, humidification hood
- Six color-coded diluter jets for different oxygen concentration:







Blue: 24%	Green: 35%
Yellow: 28%	Pink: 40%
White: 31%	Orange: 50%



Product code	Size	Secure buckle	Oxygen connector	Tubing length
GC8029C01	Adult		✓	2.1m
GC8029C02	Pediatric		✓	2.1m
GC8031C01	Adult	✓	✓	2.1m
GC8031C02	Pediatric	✓	✓	2.1m

# Tracheostomy Mask

Picture	Product code	Size	Secure buckle	Oxygen connector
	GC8029B01	Adult		✓
	GC8029B02	Pediatric		✓
	GC8029B03	Neonatal		✓
	GC8031B01	Adult	✓	✓
	GC8031B02	Pediatric	✓	✓
	GC8031B03	Neonatal	✓	✓

Picture	Product code	Size	Secure buckle	Oxygen connector
	GC8029A01	Adult		
	GC8029A02	Pediatric		
	GC8029A03	Neonatal		
	GC8031A01	Adult	✓	
	GC8031A02	Pediatric	✓	
	GC8031A03	Neonatal	✓	

# 02 AIRWAY MANAGEMENT

## Berman Airway

- Beveled tip for comfortable and gentle insertion
- Big curved airway, exactly conforming with oral structure
- Tip of pipe can reach the laryngeal correctly and make sure smooth airflow
- With or without vented septum



Vented	Non-Vented	Size
GC8023B04	GC8023D04	40mm
GC8023B05	GC8023D05	50mm
GC8023B06	GC8023D06	60mm
GC8023B07	GC8023D07	70mm
GC8023B08	GC8023D08	80mm
GC8023B09	GC8023D09	90mm
GC8023B10	GC8023D10	100mm
GC8023B11	GC8023D11	110mm



Vented	Color	Non-Vented	Color	Size
GC8023A04	Pink	GC8023C04	Pink	40mm
GC8023A05	Blue	GC8023C05	Blue	50mm
GC8023A06	Black	GC8023C06	Black	60mm
GC8023A07	White	GC8023C07	White	70mm
GC8023A08	Green	GC8023C08	Green	80mm
GC8023A09	Yellow	GC8023C09	Yellow	90mm
GC8023A10	Red	GC8023C10	Red	100mm
GC8023A11	Orange	GC8023C11	Orange	110mm

## Guedel Airway

- Big curved airway, better conforming with oral structure, safer than traditional small curve
- Color-coded insert with one piece construction, avoiding the insert pulled out by external force
- Built-in bite block



Guedel airway	Color	Size
GC8023G04	<span style="color: pink;">■</span> Pink	40mm
GC8023G05	<span style="color: blue;">■</span> Blue	50mm
GC8023G06	<span style="color: black;">■</span> Black	60mm
GC8023G07	<span style="color: white;">□</span> White	70mm
GC8023G08	<span style="color: green;">■</span> Green	80mm
GC8023G09	<span style="color: yellow;">■</span> Yellow	90mm
GC8023G10	<span style="color: red;">■</span> Red	100mm
GC8023G11	<span style="color: orange;">■</span> Orange	110mm
GC8023G12	<span style="color: purple;">■</span> Purple	120mm

## Nasopharyngeal Airway

- Soft material
- Inner striations for smooth suction catheter insertion
- Beveled tip for comfortable and gentle insertion
- Available in full sizes



Product code	Size
GC8024A12	12FR
GC8024A14	14FR
GC8024A16	16FR
GC8024A18	18FR
GC8024A20	20FR
GC8024A22	22FR
GC8024A24	24FR
GC8024A26	26FR
GC8024A28	28FR
GC8024A30	30FR
GC8024A32	32FR
GC8024A34	34FR
GC8024A36	36FR
GC8024A38	38FR

## Robertazzi Nasopharyngeal Airway

- Soft and flexible construction for patient comfort
- Round tip allows for gentle insertion



Product code	Size
GC8024C12	12FR
GC8024C14	14FR
GC8024C16	16FR
GC8024C18	18FR
GC8024C20	20FR
GC8024C22	22FR
GC8024C24	24FR
GC8024C26	26FR
GC8024C28	28FR
GC8024C30	30FR
GC8024C32	32FR
GC8024C34	34FR
GC8024C36	36FR
GC8024C38	38FR

# Endotracheal Tube

- All tubes are fitted with 15mm standard connector, which insures compatibility with circuit connectors
- Smooth tip aids atraumatic intubations
- Self-closing valve with Luer-lock syringe attachment piece
- Murphy eye



Uncuffed	Cuffed	Size
GC8039A020	GC8039B020	2.0
GC8039A025	GC8039B025	2.5
GC8039A030	GC8039B030	3.0
GC8039A035	GC8039B035	3.5
GC8039A040	GC8039B040	4.0
GC8039A045	GC8039B045	4.5
GC8039A050	GC8039B050	5.0
GC8039A055	GC8039B055	5.5
GC8039A060	GC8039B060	6.0
GC8039A065	GC8039B065	6.5
GC8039A070	GC8039B070	7.0
GC8039A075	GC8039B075	7.5
GC8039A080	GC8039B080	8.0
GC8039A085	GC8039B085	8.5
GC8039A090	GC8039B090	9.0
GC8039A095	GC8039B095	9.5
GC8039A100	GC8039B100	10.0

# Reinforced Endotracheal Tube

- Soft, flexible PVC tube with spiral-wound reinforcing wire is designed to help reduce risk of kinking
- Reinforcing wire is sealed tightly against bonded connector
- Full magill curve
- Murphy eye
- Latex free



Uncuffed	Cuffed	Size
GC8040A030	GC8040B030	3.0
GC8040A035	GC8040B035	3.5
GC8040A040	GC8040B040	4.0
GC8040A045	GC8040B045	4.5
GC8040A050	GC8040B050	5.0
GC8040A055	GC8040B055	5.5
GC8040A060	GC8040B060	6.0
GC8040A065	GC8040B065	6.5
GC8040A070	GC8040B070	7.0
GC8040A075	GC8040B075	7.5
GC8040A080	GC8040B080	8.0
GC8040A085	GC8040B085	8.5
GC8040A090	GC8040B090	9.0
GC8040A095	GC8040B095	9.5

# Suction Plus Endotracheal Tube


- For airway management by oral/nasal intubation of the trachea and for evacuation or drainage of secretion from the subglottic space
- Helps decreasing the rate of ventilator associated pneumonia(VAP) in the hospital and to reduce VAP related costs
- Integral suction lumen allows continuous suctioning without risking trauma to the vocal cords as with manual catheter suctioning



Standard suction plus	Reinforced suction plus	Size
GC8039C060	GC8040C060	6.0
GC8039C065	GC8040C065	6.5
GC8039C070	GC8040C070	7.0
GC8039C075	GC8040C075	7.5
GC8039C080	GC8040C080	8.0
GC8039C085	GC8040C085	8.5
GC8039C090	GC8040C090	9.0

# Oral Preformed Endotracheal Tube

- Designed with preformed curvature to reduce tube kinking
- Preformed shape directs tube downward, to rest on patients chin
- Radiopaque x-ray line
- Murphy eye



Uncuffed	Cuffed	Size
GC8042A030	GC8042B030	3.0
GC8042A035	GC8042B035	3.5
GC8042A040	GC8042B040	4.0
GC8042A045	GC8042B045	4.5
GC8042A050	GC8042B050	5.0
GC8042A055	GC8042B055	5.5
GC8042A060	GC8042B060	6.0
GC8042A065	GC8042B065	6.5
GC8042A070	GC8042B070	7.0
GC8042A075	GC8042B075	7.5
GC8042A080	GC8042B080	8.0
GC8042A085	GC8042B085	8.5
GC8042A090	GC8042B090	9.0
GC8042A095	GC8042B095	9.5
GC8042A100	GC8042B100	10.0

# Nasal Preformed Endotracheal Tube

- Preformed curvature reduce tube kinking
- The preformed shape directs tube over patients forehead, reducing pressure on nares
- Radiopaque x-ray line
- Murphy eye



Uncuffed	Cuffed	Size
GC8043A030	GC8043B030	3.0
GC8043A035	GC8043B035	3.5
GC8043A040	GC8043B040	4.0
GC8043A045	GC8043B045	4.5
GC8043A050	GC8043B050	5.0
GC8043A055	GC8043B055	5.5
GC8043A060	GC8043B060	6.0
GC8043A065	GC8043B065	6.5
GC8043A070	GC8043B070	7.0
GC8043A075	GC8043B075	7.5
GC8043A080	GC8043B080	8.0
GC8043A085	GC8043B085	8.5
GC8043A090	GC8043B090	9.0
GC8043A095	GC8043B095	9.5
GC8043A100	GC8043B100	10.0




## Right/Left-sided Endobronchial Tube

- Endobronchial tube, connector set and two control valve
- Suction catheters with graduation mark
- Made of thermosensitive, clear PVC
- Two soft low-pressure cuffs
- Predesigned stylet
- Softer, straighter tube

Left-sided	Right-sided	Size
GC8051A028	GC8051B028	28FR
GC8051A032	GC8051B032	32FR
GC8051A035	GC8051B035	35FR
GC8051A037	GC8051B037	37FR
GC8051A039	GC8051B039	39FR
GC8051A041	GC8051B041	41FR

## CELESTIAL Adjustable Tracheostomy Tube

- The tube body is printed with markings to help medical staff position the tube accurately
- The protective wings are adjustable, soft, and smooth, making it comfortable for the patient to wear
- Supplied with cotton neck tape



Uncuffed	Cuffed	Tube size	Tube I.D.	Tube O.D.	Tube length
GC8041C060N	GC8041D060N	6.0mm	6.0mm	8.3mm	99.0mm
GC8041C070N	GC8041D070N	7.0mm	7.0mm	9.6mm	102.0mm
GC8041C080N	GC8041D080N	8.0mm	8.0mm	11.0mm	112.5mm
GC8041C090N	GC8041D090N	9.0mm	9.0mm	12.4mm	118.5mm
GC8041C100N	GC8041D100N	10.0mm	10.0mm	13.8mm	131.5mm

## CELESTIAL Tracheostomy Tube with NeoSuction

The NeoSuction tracheostomy has an additional subglottic port which has 2 possible functions:

- Removes pooled subglottic secretions using above-cuff suctioning
- Enables above-cuff voicing (if patient has been assessed as suitable by a Speech Pathologist).



Cuffed	Tube size	Tube I.D.	Tube O.D.	Tube length
GC8041E060N	6.0mm	6.0mm	8.3mm	99.0mm
GC8041E070N	7.0mm	7.0mm	9.6mm	102.0mm
GC8041E080N	8.0mm	8.0mm	11mm	112.5mm
GC8041E090N	9.0mm	9.0mm	12.4mm	118.5mm
GC8041E100N	10.0mm	10.0mm	13.8mm	131.5mm

# CELESTIAL Reinforced Tracheostomy Tube



Cuffed	Tube size	Tube I.D.	Tube O.D.	Tube length
GC8041F060N	6.0mm	6.0mm	8.3mm	99.0mm
GC8041F070N	7.0mm	7.0mm	9.6mm	102.0mm
GC8041F080N	8.0mm	8.0mm	11.0mm	112.5mm
GC8041F090N	9.0mm	9.0mm	12.4mm	118.5mm
GC8041F100N	10.0mm	10.0mm	13.8mm	131.5mm

# SERENITY Tracheostomy Tube



Uncuffed	Cuffed	Tube size	Tube I.D.	Tube O.D.	Tube length
GC8041R060N	GC8041S060N	6.0mm	6.0mm	9.2mm	64.5mm
GC8041R065N	GC8041S065N	6.5mm	6.5mm	9.8mm	67.0mm
GC8041R070N	GC8041S070N	7.0mm	7.0mm	10.5mm	70.0mm
GC8041R075N	GC8041S075N	7.5mm	7.5mm	11.3mm	73.0mm
GC8041R080N	GC8041S080N	8.0mm	8.0mm	11.9mm	75.5mm
GC8041R085N	GC8041S085N	8.5mm	8.5mm	12.6mm	78.0mm
GC8041R090N	GC8041S090N	9.0mm	9.0mm	13.3mm	81.0mm
GC8041R100N	GC8041S100N	10.0mm	10.0mm	14.0mm	87.5mm

# SERENITY Tracheostomy Tube Kit




Uncuffed	Cuffed	Tube size	Tube I.D.	Tube O.D.	Tube length
GC8041R060C	GC8041S060C	6.0mm	6.0mm	9.2mm	64.5mm
GC8041R065C	GC8041S065C	6.5mm	6.5mm	9.8mm	67.0mm
GC8041R070C	GC8041S070C	7.0mm	7.0mm	10.5mm	70.0mm
GC8041R075C	GC8041S075C	7.5mm	7.5mm	11.3mm	73.0mm
GC8041R080C	GC8041S080C	8.0mm	8.0mm	11.9mm	75.5mm
GC8041R085C	GC8041S085C	8.5mm	8.5mm	12.6mm	78.0mm
GC8041R090C	GC8041S090C	9.0mm	9.0mm	13.3mm	81.0mm
GC8041R100C	GC8041S100C	10.0mm	10.0mm	14.0mm	87.5mm

## SERENITY Tracheostomy Tube with NeoSuction



Product code	Tube size	Tube I.D.	Tube O.D.	Tube length
GC8041T060N	6.0mm	6.0mm	9.2mm	64.5mm
GC8041T070N	7.0mm	7.0mm	10.5mm	70.0mm
GC8041T075N	7.5mm	7.5mm	11.3mm	73.0mm
GC8041T080N	8.0mm	8.0mm	11.9mm	75.5mm
GC8041T085N	8.5mm	8.5mm	12.6mm	78.0mm
GC8041T090N	9.0mm	9.0mm	13.3mm	81.0mm
GC8041T100N	10.0mm	10.0mm	14.0mm	87.5mm

## SERENITY Tracheostomy Tube Kit with NeoSuction



Product code	Tube size	Tube I.D.	Tube O.D.	Tube length
GC8041T060C	6.0mm	6.0mm	9.2mm	64.5mm
GC8041T070C	7.0mm	7.0mm	10.5mm	70.0mm
GC8041T075C	7.5mm	7.5mm	11.3mm	73.0mm
GC8041T080C	8.0mm	8.0mm	11.9mm	75.5mm
GC8041T085C	8.5mm	8.5mm	12.6mm	78.0mm
GC8041T090C	9.0mm	9.0mm	13.3mm	81.0mm
GC8041T100C	10.0mm	10.0mm	14.0mm	87.5mm

## GRASMERE Fenestrated Tracheostomy Tube



Uncuffed	Cuffed	Tube I.D.	Tube O.D.	Tube length
GC8041G060N	GC8041H060N	6.0mm	9.2mm	64.5mm
GC8041G070N	GC8041H070N	7.0mm	10.5mm	70.0mm
GC8041G075N	GC8041H075N	7.5mm	11.3mm	73.0mm
GC8041G080N	GC8041H080N	8.0mm	11.9mm	75.5mm
GC8041G085N	GC8041H085N	8.5mm	12.6mm	78.0mm
GC8041G090N	GC8041H090N	9.0mm	13.3mm	81.0mm
GC8041G100N	GC8041H100N	10.0mm	14.0mm	87.5mm

## GRASMERE Fenestrated Tracheostomy Tube Kit



Uncuffed	Cuffed	Tube size	Tube I.D.	Tube O.D.	Tube length
GC8041G060B	GC8041H060B	6.0mm	6.0mm	9.2mm	64.5mm
GC8041G070B	GC8041H070B	7.0mm	7.0mm	10.5mm	70.0mm
GC8041G075B	GC8041H075B	7.5mm	7.5mm	11.3mm	73.0mm
GC8041G080B	GC8041H080B	8.0mm	8.0mm	11.9mm	75.5mm
GC8041G085B	GC8041H085B	8.5mm	8.5mm	12.6mm	78.0mm
GC8041G090B	GC8041H090B	9.0mm	9.0mm	13.3mm	81.0mm
GC8041G100B	GC8041H100B	10.0mm	10.0mm	14.0mm	87.5mm

## AZURE Non-adjustable Tracheostomy Tube



Uncuffed	Cuffed	Size
GC8041A060N	GC8041B060N	6.0
GC8041A065N	GC8041B065N	6.5
GC8041A070N	GC8041B070N	7.0
GC8041A075N	GC8041B075N	7.5
GC8041A080N	GC8041B080N	8.0
GC8041A085N	GC8041B085N	8.5
GC8041A090N	GC8041B090N	9.0
GC8041A095N	GC8041B095N	9.5
GC8041A100N	GC8041B100N	10.0



## Intubating Stylet

- Made of malleable aluminum covered with plastic sheath, which has a lubricated Satin-Slip surface
- Helps reduce friction between stylet and tracheal tube to allow for easy insertion and withdrawal
- Sheath extends beyond tip to reduce the risk of trauma



Product code	Size
GC8048A06	06FR
GC8048A10	10FR
GC8048A12	12FR
GC8048A14	14FR

# Endotracheal Tube Introducer

- Bougie type introducer
- Flexible reinforcement
- Graduation mark at each side
- Angled tip for ease of intubation
- Available with solid type and hollow type




Hollow type	DEHP free hollow type	Size
GC8052A06	/	06FR
GC8052A10	GC8052B10	10FR
GC8052A15	GC8052B15	15FR

# Disposable Laryngeal Mask Airway



PVC/Straight open	Silicone/Straight open	Silicone/Glottic bars	Size	Patient weight
GC8001A10	GC8001B10	GC8001G10	1	<5kg
GC8001A15	GC8001B15	GC8001G15	1.5	5-10kg
GC8001A20	GC8001B20	GC8001G20	2	10-20kg
GC8001A25	GC8001B25	GC8001G25	2.5	20-30kg
GC8001A30	GC8001B30	GC8001G30	3	30-50kg
GC8001A40	GC8001B40	GC8001G40	4	50-70kg
GC8001A50	GC8001B50	GC8001G50	5	>70kg

# 90° Preformed Laryngeal Mask Airway



Product code	Size	Patient weight	Open type
GC8001K10	1	<5kg	Straight open
GC8001K15	1.5	5-10kg	Straight open
GC8001K20	2	10-20kg	Straight open
GC8001K25	2.5	20-30kg	Straight open
GC8001K30	3	30-50kg	Straight open
GC8001K40	4	50-70kg	Straight open
GC8001K50	5	>70kg	Straight open

# Reinforced Laryngeal Mask Airway



Silicone/Straight open	Silicone/Glottic bars	PVC/Straight open	Size	Patient weight
GC8001D10	GC8001H10	GC8001C10	1	<5kg
GC8001D15	GC8001H15	GC8001C15	1.5	5-10kg
GC8001D20	GC8001H20	GC8001C20	2	10-20kg
GC8001D25	GC8001H25	GC8001C25	2.5	20-30kg
GC8001D30	GC8001H30	GC8001C30	3	30-50kg
GC8001D40	GC8001H40	GC8001C40	4	50-70kg
GC8001D50	GC8001H50	GC8001C50	5	>70kg


# Reusable Laryngeal Mask Airway



Silicone/Straight open	Reinforced reusable	Size	Patient weight	Open type
GC8001E10	GC8001F10	1	<5kg	Glottic bars
GC8001E15	GC8001F15	1.5	5-10kg	Glottic bars
GC8001E20	GC8001F20	2	10-20kg	Glottic bars
GC8001E25	GC8001F25	2.5	20-30kg	Glottic bars
GC8001E30	GC8001F30	3	30-50kg	Glottic bars
GC8001E40	GC8001F40	4	50-70kg	Glottic bars
GC8001E50	GC8001F50	5	>70kg	Glottic bars

# Double-lumen Silicone Laryngeal Mask Airway


- For establishing an artificial airway in emergency situations
- Silicone, double lumen, for single use



Product code	Size	Patient weight
GC8001S10	1	<5kg
GC8001S15	1.5	5-10kg
GC8001S20	2	10-20kg
GC8001S25	2.5	20-30kg
GC8001S30	3	30-50kg
GC8001S40	4	50-70kg
GC8001S50	5	>70kg

# 03 RESPIRATORY CARE

## Mucus Trap with Vacuum Breaker



Product code	Description
GC8054C08	08FR, 20ml
GC8054C10	10FR, 20ml

# Mucus Trap with Mouthpiece



Product code	Description
GC8054K07	07FR, 20ml
GC8054K08	08FR, 20ml
GC8054K10	10FR, 20ml

# Mucus Trap with Tapered Connector



Product code

GC8054L14

Description

14FR, 20ml


# Mucus Extractor

- Transparent chamber for visual examination of the content
- Catheter is equipped with a round, open tip for atraumatic intubation
- The spare cap is provided to seal the container for safe transportation of specimen to the laboratory or aseptic disposal of the container



GC8054D00

Product code	Description
GC8054D06	06FR, 25ml, open tip
GC8054D08	08FR, 25ml, open tip
GC8054D10	10FR, 25ml, open tip
GC8054D12	12FR, 25ml, open tip
GC8054D14	14FR, 25ml, open tip
GC8054D00	14FR, 25ml, funnel connector



**GC8054E00**

Product code	Description
GC8054E10	10FR, 25ml, without side eye
GC8054E12	12FR, 25ml, without side eye
GC8054E14	14FR, 25ml, without side eye
GC8054E00	14FR, 25ml, with gradual step connector for funnel suction catheters

# Mucus Specimen Trap

- Mucus specimen traps allow for easy, convenient collection of mucus specimen samples
- After specimen collection, the two connectors are joined to close the system



Product code	Description
GC8057A04	40ml, with spare screw cap, removal adaptor & patient label
GC8057A07	70ml, with spare screw cap, removal adaptor & patient label



Product code	Description
GC8057B04	40ml, with spare screw cap & patient label
GC8057B07	70ml, with spare screw cap & patient label

# Mucus Specimen Traps

Designed for medical convenience, it features a transparent container for clear visibility, a leak-proof design for safety, and easy-to-use connectors for seamless integration with suction devices. Ideal for respiratory and sputum samples, it ensures precise specimen collection for accurate diagnostic results.



Product code	Description
GC8057C07	70ml, with 11cm tubing, for use in bronchoscopy
GC8054M00C01	25ml, with reduced length connection to suction catheter to ease movement of secretions to the trap
GC8054M00C02	25ml, with 11cm tubing, for use in bronchoscopy
GC8054N00C01	25ml, for use in bronchoscopy

## Suction Catheter


GCMedica has a various range of suction catheter. Genuine factory allows us to manufacture a range of lengths, connector types, sizes and tips to meet all clinical requirements.

We have standard tip, milly tip, autramatic tip and curved tip.

Connector are available with Funnel connector and Transparent Y connector with vacuum control.



STANDARD TIP



MILLY TIP



ATRAUMATIC TIP



CURVED TIP

## Suction Catheter - Transparent Y Connector with Vacuum Control


Size	06FR	08FR	10FR	12FR	14FR	16FR	18FR
<b>Standard Tip</b>	GC8018Y06	GC8018Y08	GC8018Y10	GC8018Y12	GC8018Y14	GC8018Y16	GC8018Y18
<b>Milly Tip</b>	GC8018F06	GC8018F08	GC8018F10	GC8018F12	GC8018F14	GC8018F16	GC8018F18
<b>Autramatic Tip</b>	GC8018M06	GC8018M08	GC8018M10	GC8018M12	GC8018M14	GC8018M16	GC8018M18
<b>Curved Tip</b>	GC8018B06	GC8018B08	GC8018B10	GC8018B12	GC8018B14	GC8018B16	GC8018B18

## Suction Catheter - Funnel Connector

Size	06FR	08FR	10FR	12FR	14FR	16FR	18FR
<b>Standard Tip</b>	GC8018P06	GC8018P08	GC8018P10	GC8018P12	GC8018P14	GC8018P16	GC8018P18
<b>Milly Tip</b>	GC8018I06	GC8018I08	GC8018I10	GC8018I12	GC8018I14	GC8018I16	GC8018I18
<b>Autramatic Tip</b>	GC8018L06	GC8018L08	GC8018L10	GC8018L12	GC8018L14	GC8018L16	GC8018L18
<b>Curved Tip</b>	GC8018K06	GC8018K08	GC8018K10	GC8018K12	GC8018K14	GC8018K16	GC8018K18

# MDI Spacer

Assistive device for inhalation of aerosolized medication for asthma treatment.



Product code	Description
GC8064A01	300ml, Built-in speaker

# 04 ANESTHESIA



Double Swivel Wye  
Connector 360° Shoulder



## Corrugated Breathing Circuit

### Adult

Product code	Wye connector	Length	BV filter	HME	Gas sampling line	Anesthesia mask	Breathing bag
GC40COR22	Parallel wye	40"					
GC60COR22	Parallel wye	60"					
GC72COR22	Parallel wye	72"					
GC40CORGSP22	Parallel wye with port	40"					
GC60CORGSP22	Parallel wye with port	60"					
GC72CORGSP22	Parallel wye with port	72"					
GC40CORSW22	Swivel parallel wye with port	40"					
GC60CORSW22	Swivel parallel wye with port	60"					
GC72CORSW22	Swivel parallel wye with port	72"					
GC60CORBW22	Bifurcated wye	60"					
GC80CA1A01	Parallel wye	40"	✓				
GC80CA1B01	Parallel wye	40"		✓	✓		3L
GC80CA1C01	Parallel wye	40"	2				3L
GC80CA1D01	Parallel wye	40"					0.5L
GC80CA1E01	Parallel wye	40"				5#	2L
GC80CB2A01	Parallel wye with port	60"	✓				
GC80CC4B02	Parallel wye with port	90"	✓		✓		3L
GC80CB2C01	Swivel Parallel wye with port	60"	2				3L
GC80CA2E01	Parallel wye	60"				5#	2L

### Pediatric

Product code	Wye connector	Length	BV filter	HME	Gas sampling line	Anesthesia mask	Breathing bag
GC60COR15	Parallel wye	60"					
GC72COR15	Parallel wye	72"					
GC60CORGSP15	Parallel wye with port	60"					
GC72CORGSP15	Parallel wye with port	72"					
GC80CF2D01	Parallel wye with port	60"					0.5L

# Expandable Breathing Circuit



## Adult

Product code	Wye connector	Length	BV filter	HME	Gas sampling line	Anesthesia mask	Breathing bag
GC60EXP22	Parallel wye	60"					
GC72EXP22	Parallel wye	72"					
GC90EXP22	Parallel wye	90"					
GC90EXPGSP22	Parallel wye with port	90"					
GC90EXPSW22	Swivel parallel wye with port	90"					
GC90EXPBW22	Bifurcated wye with port	90"					
GC80EB4F02	Parallel wye	90"		√		5#	3L
GC80EC4F03	Swivel parallel wye with port	90"	√			5#	3L
GC80EC3G01	Swivel parallel wye with port	72"	√		√	5#	3L
GC80EA3G03	Parallel wye	72"	2		√	5#	3L
GC80EA1D02	Parallel wye	40"					1L

## Pediatric

Product code	Wye connector	Length	BV filter	HME	Gas sampling line	Anesthesia mask	Breathing bag
GC60EXP15	Parallel wye	60"					
GC90EXP22	Parallel wye	90"					
GC108EXP22	Parallel wye	108"					
GC120EXP22	Parallel wye	120"					
GC60EXPGSP15	Parallel wye with port	60"					
GC90EXPGSP15	Parallel wye with port	90"					

# Coaxial Breathing Circuit

- Inner and outer tube design. Attaches via standard connections
- Circuit tubing are available with adult and pediatric corrugated circuit
- Wye connector are available with parallel wye with or without gas sampling port, swivel wye, bifurcated wye connector



Adult/Tube I.D.22mm	Pediatric/Tube I.D.15mm	Length
GC40COAX	/	40"
GC60COAX	GC60COAX15	60"
GC72COAX	GC72COAX15	72"
GC90COAX	/	90"
GC120COAX	/	120"

# Smoothbore Breathing Circuit



Product code	Length	Tube I.D.
GC851101	1.6m	22mm
GC851102	1.6m	15mm

# Catheter Mount

- With highly lightweight flexible tube, which prevents lateral pressure and minimizes drag from breathing circuit
- Available with double swivel elbow for more convenience while operating



Product code	Adapter	Tube length	Tube type
GC8036B01	Double swivel elbow	150mm	Expandable tube
GC8036B02	Double swivel elbow	150mm	Corrugated tube
GC8036A01	Standard elbow	150mm	Expandable tube
GC8036A02	Standard elbow	150mm	Corrugated tube
GC8036C02	Straight connector	150mm	Flexible tube with HME

## Anesthesia Mask

- Unique ergonomic shape fits all facial configurations for improved sealing capabilities
- Soft, pliable air filled cushion with and without an inflation valve

Product code	Hook color	Connection port	Intend for	Size
GC8020A01	White	15mm OD	Neonatal	0#
GC8020A02	Pink	15mm OD	Infant	1#
GC8020A03	Yellow	22mm OD	Pediatric	2#
GC8020A04	Green	22mm OD	Adult-S	3#
GC8020A05	Red	22mm OD	Adult-M	4#
GC8020A06	Blue	22mm OD	Adult-L	5#
GC8020A07	Orange	22mm OD	Adult-XL	6#



## Anes-Air Cushion Mask

- Complying with biocompatibility standards, soft air cushion ensures flexibility, air tightness and low odor
- The mask is specially designed to fit comfortable on patients face
- Inflating air cushion for better sealing. High transparency allows good visibility

Product code	Hook color	Connection port	Intend for	Size
GC8020B01	/	15mm OD	Neonatal	0#
GC8020B02	/	15mm OD	Infant	1#
GC8020B03	/	15mm OD	Pediatric	2#
GC8020B04	Green	22mm OD	Adult-S	3#
GC8020B05	Red	22mm OD	Adult-M	4#
GC8020B06	Blue	22mm OD	Adult-L	5#
GC8020B07	Orange	22mm OD	Adult-XL	6#



## Eco-Anes Mask

- Ergonomically designed with easy use grip for practicality
- High transparency shell allows good visibility
- Better sealing and light weight without the need to inflate the cushion

Product code	Cushion color	Connection port	Intend for	Size
GC8020D01	Blue	15mm OD	Infant	1#
GC8020D02	Orange	15mm OD	Pediatric	2#
GC8020D03	Pink	22mm OD	Adult-S	3#
GC8020D04	White	22mm OD	Adult-M	4#
GC8020D05	Yellow	22mm OD	Adult-L	5#
GC8020D06	Green	22mm OD	Adult-XL	6#



# Manual Resuscitator

Available with resuscitator kit, including manual resuscitator bag, anesthesia mask, oral airway, oxygen tube



PVC manual resuscitator	Silicone manual resuscitator	Size	Kit
GC8044A01	GC8044B01	Adult	Manual resuscitator only
GC8044A02	GC8044B02	Pediatric	Manual resuscitator only
GC8044A03	GC8044B03	Neonatal	Manual resuscitator only
GC8045A01	GC8045B01	Adult	Resuscitator kit
GC8045A02	GC8045B02	Pediatric	Resuscitator kit
GC8045A03	GC8045B03	Neonatal	Resuscitator kit

# Incentive Breathing Exerciser

- It helps the patient to recover normal respiration after a chest or abdominal surgery
- It encourages deep and prolonged inspiration
- Three chamber design allows the patient to lift the balls without any resistance to achieve the maximum volume in minimum time
- Available with adjustable and non-adjustable type



Product code	Volume	Size
GC8002A01	600-1200cc/second	Non-adjustable
GC8002B01	600-1200cc/second	Adjustable

# Volumetric Breathing Exerciser

- Two sizes available: 5000ml and 2500ml
- Colorful volume labeling facilitating patient evaluation and use
- Clean flow window for convenient observation




Product code	Type	Volume
GC8002C03	Non-Adjustable	Up to 5000cc
GC8002C02	Non-Adjustable	Up to 2500cc

## Gas Sampling Line



Product code	Length	Type	Volume
GC800301	3m	1.2mm	Male + Male
GC800302	3m	1.2mm	Male + Female

# Heat Moisture Exchange (HME)




High efficient HME


Product code	GC800702
Moisture output	31mg/H <sub>2</sub> O
Luer lock port	√
Compressible volume	45ml
Tide volume	150-1500ml
Resistance(cmH <sub>2</sub> O)60LPM	60pa




Eco-mini HME



GC800705




GC800709



GC800713

Product code	GC800705	GC800707	GC800709	GC800711	GC800713	GC800715
Moisture output	30mg/H <sub>2</sub> O	31mg/H <sub>2</sub> O	30mg/H <sub>2</sub> O	31mg/H <sub>2</sub> O	30mg/H <sub>2</sub> O	31mg/H <sub>2</sub> O
Luer lock port	√	√	√	√	√	√
Compressible volume	40ml	30ml	40ml	30ml	40ml	30ml
Tide volume	150-1500ml	900-1500ml	150-1500ml	900-1500ml	150-1500ml	900-1500ml
Resistance(cmH <sub>2</sub> O)30LPM		50pa		50pa		50pa
Resistance(cmH <sub>2</sub> O)60LPM	60pa		60pa		250pa	
Angled					√	
Size	Adult	Pediatric	Adult	Pediatric	Adult	Pediatric

# Heat Moisture Exchange Filter (HMEF)




High efficient HMEF


Product code	GC800703
Moisture output	31mg/H <sub>2</sub> O
Luer lock port	√
Compressible volume	45ml
Tide volume	150-1500ml
Filtration efficiency	>99.99%
Resistance(cmH <sub>2</sub> O)60LPM	240pa




Eco-mini HMEF



GC800705




GC800709



GC800713

Product code	GC800706	GC800708	GC800710	GC800712	GC800714	GC800716
Moisture output	30mg/H <sub>2</sub> O	31mg/H <sub>2</sub> O	30mg/H <sub>2</sub> O	31mg/H <sub>2</sub> O	30mg/H <sub>2</sub> O	31mg/H <sub>2</sub> O
Luer lock port	√	√			√	√
Compressible volume	40ml	30ml	40ml	30ml	40ml	30ml
Tide volume	150-1500ml	900-1500ml	150-1500ml	900-1500ml	150-1500ml	900-1500ml
Filtration efficiency	>99.99%	>99.99%	>99.99%	>99.99%	>99.99%	>99.99%
Resistance(cmH <sub>2</sub> O)30LPM	60pa		255pa		66pa	
Angled					√	
Size	Adult	Pediatric	Adult	Pediatric	Adult	Pediatric

# Bacterial/Viral Filter



<b>Product code</b>	<b>GC800701</b>
Luer lock port	√
Tide volume	150-1500ml
Filtration efficiency	>99.99%
Resistance(cmH <sub>2</sub> O)60LPM	230pa

# 05 ABC MOUTHPIECE AND FILTER KIT

## ABC Mouthpiece and Filter Kit

The mouthpiece and filter kit is intended to be used with active breathing coordinator for DIBH (Deep Inspiration Breathing Hold), a technique to treat cancer in the breast, chest wall or occasionally the abdomen. The mouth and nose clip help patient inhale and exhale via mouth and take a deep breath during radiation treatment under doctor's guide.



Product code	Description
GC8049A01	A kit including mouth piece, 17" breathing hose, filter and nose clip.



## CONTACT US

---



**LOUJIN INDUSTRIAL PARK, SHUOFANG,  
WUXI 214143, Jiangsu, P.R.CHINA**



**86-510-82733575**



**INFO@GCMEDICA.COM**



**WWW.GCMEDICA.COM**