



UROLOGY

GCMEDICA ENTERPRISE LTD.(WUXI)
GCMEDICA ENTERPRISE LTD.(WUXI)

01

URINARY DEVICE

√ UROSTREAM DRAIN BAG	01
√ URINE DRAINAGE BAG	02
Basic Urine Bag	02
Low-Bed Urine Bag	02
Urine Meters	03-05
Balsam Urine Bag	06
Pear-Shaped Drainage Bag	07
√ LEG BAG	08
√ TRIANGULAR GRADUATES	08
√ MALE URINAL	09
√ UNISEX URINAL	09

UROSTREAM DRAIN BAG

Urology Drain Bag is designed to fit most hoop frames. This unique design includes a vinyl flap that goes under the patient to direct fluids into the bag, a flexible 6' long hose, and a tissue filter to allow for the collection of the tissue specimens. This optional design makes things safer and easier for your OR personnel by allowing them to utilize your existing suction equipment. The result is better management and control of infectious fluid waste.

FEATURES:


- Mounts with integral clamps to the side rail
- A flexible band provides excellent access while providing superior fluid containment
- Compatible with all suction systems for exact liquid collection
- Open mesh filter allow for fast drainage of fluids
- Durable, drape has flexible hose
- Velcro securement eliminates the need to use tape



Product code	Description
GC05701	Fixed filter
GC05702	Removable filter

URINE DRAINAGE BAG

Basic Urine Bag



Push & pull valve



Twist valve

Product code	Description
GC051A00	2000ml, no drain valve, anti-reflux valve
GC051A01	2000ml, push drain valve, anti-reflux valve
GC051A03	2000ml, twist drain valve, anti-reflux valve
GC051A12	2000ml, slide drain valve, anti-reflux valve
GC051A22	2000ml, slide drain valve, without anti-reflux valve

Low-Bed Urine Bag

Low-Bed Drainage Bag features to facilitate proper drainage and anti-reflux valve for back-flow prevention. The exclusive slide-tap drainage port allows for minimal splashing during drainage. The drainage bag has a reversible design for versatile placement on either side of the bed.

Product code GC051M01

FEATURES:

- With cross valve
- With anti-reflux valve
- With reinforced double hanger & rope hanger



Urine Meters

Drainage Bag with Urine Meter is a closed suction system used to accurately measure urine output and collect the urine. The device features a rigid, clear and transparent measured volume chamber for critically ill and surgical patients with 1ml graduations for easy reading of output measurements.

FEATURES:

- Needleless sample port with blunt cannula access
- Adjustable bed sheet clamp for securing drain bag tubing
- Designed with anti-reflux valve



GC051E11



GC051G01




GC051G02

Product code	Description
GC051G01	400ml Urine Meter with 2,500ml Drain Bag, 150cm tubing
GC051G02	400ml Urine Meter with 2,500ml Drain Bag, Metal Clamp, Slide Tap, 150cm tubing
GC051E11	500ml Urine Meter with 2,600ml Drain Bag, Slide Tap, 150cm tubing

FEATURES:


- Detachable urinary bag with urine meter on top and urine bag on bottom
- Rotatable meter cap that controls the urine switch
- Reinforced double hanger, rope hanger and double transparent straps
- Special polygonal urine meter
- Flip-type outlet valve



Product code	Description
GC051J01	500ml Urine Meter with 2000ml Drain bag, Rotating Tap, 120cm tubing

FEATURES:

- Dual No-Return Valves:
Minimizes the risk of retrograde infection for safer use
- Double Lumen Tube:
Reduces urine adhesion to the tube, ensuring free and uninterrupted urine flow
- High-Performance Bacterial Filter
- Easy Collection Bag Replacement
- Ergonomic Handle for Easy Connection




Product code	Tube size	Volume	Description
GC051P01NFC01	150cm	500ml metered	Needle-free sample port (Incision)
GC051P01NFU01	150cm	chamber + 2000ml	Needle-free sample port (Incision), PVC Free
GC051P02NFC01	150cm	bag	Needle-free sample port (Luer)

Balsam Urine Bag

FEATURES:

- 2000ml volume
- Diamond shaped bag
- Inlet anti-reflux valve
- With reinforced hanger
- With bed sheet clamp
- Φ 10.0x110cm inlet tubing
- Needleless sample port with cap
- Slide T-tap outlet valve
- Latex-free, single packed



Product code	Description
GC051B01	2000ml, with 110cm tube, slide tap valve

Pear-Shaped Drainage Bag

FEATURES:

- Pear-Shaped bag
- With reinforced single hanger & rope hanger



Product code	Description
GC051C01	2000ml, φ 10.0 × 1100mm, with reinforced single hanger & rope hanger

LEG BAG

Urinary Leg Bag is durably constructed for odor containment and leak resistance. The anti-reflux chamber reduces gross reflux.

Product code	Spec
GC055B05	500ml
GC055B07	750ml
GC055B10	1000ml

FEATURES:

- Elliptic shape, with anti-reflux valve
- Design allows for even filling and draining
- Includes two flexible straps
- Single use
- Latex free



TRIANGULAR GRADUATES

FEATURES:

- Ideal for measuring intake or output, these easy-to-handle graduated containers have a frosted panel for writing pertinent patient information
- Printed black graduations for increased visibility and ease of measuring fluid level
- Graduated in 50cc and 2oz increments up to 1000cc and 32oz



Product code	Description
GC1922A01	Triangular Graduate 1000ml

MALE URINAL

The male urinal is designed to prevent spills and odors with an attached snap-on lid. It features a sturdy grip for easy handling and can be used in several positions.



Product code	Description
GC065D10	Male, 1000ml

UNISEX URINAL

The unisex urinal is both convenient and portable, and is complete with an easy carry handle, designed specifically for either male or female use via the adaptor provided. After use the urinal can be sealed using the attached anti-spill lid for safe transportation and disposal.



Product code	Description
GC065E10	Unisex, 1000ml

02


URINARY CATHETER

✓ PORTABLE INTERMITTENT CATHETER	10
CathVantage™ Portable Hydrophilic Intermittent Catheter	10
CathComfort™ Portable Lubricant Gel Intermittent Catheter	10
✓ NELATON CATHETER	11
✓ HYDROPHILIC CATHETER	12
✓ CATHETER SPIGOT	12
✓ RECTAL CATHETER	13
✓ SILICONE CATHETER	14
✓ SILICONE MALE CATHETER	15


PORTABLE INTERMITTENT CATHETER

CathVantage™ Portable Hydrophilic Intermittent Catheter

The CathVantage™ Portable Hydrophilic Intermittent Catheter is a compact and portable option that's about half the size of standard catheters and has a non-medical appearance. Its design makes it easy to store, carry, and dispose of without sacrificing usability. The catheter's coating provides even lower friction than traditional hydrophilic coatings, ensuring greater comfort and minimizing the risk of urethral micro-trauma.



Patent pending




Product code	Description	Size	Effective length
GC05808F08J		08FR	7-9CM
GC05808F10J		10FR	
GC05808F12J		12FR	
GC05808F14J		14FR	
GC05818F08J	Female	08FR	7-9CM
GC05818F10J		10FR	
GC05818F12J		12FR	
GC05818F14J		14FR	
GC05818F08K		08FR	10-15CM
GC05818F10K		10FR	
GC05818F12K		12FR	
GC05818F14K		14FR	
GC05818M12L	Male	12FR	30-40CM



CathComfort™ Portable Lubricant Gel Intermittent Catheter

FEATURES:

- Smooth eyelets for comfort
- Compact and portable design with a twist-off cap, universal funnel for broad compatibility
- Catheter Kits are available




Pre-lubricated



Product code	Description	Size	Effective length
GC05828F08J	Female	08FR	7-9CM
GC05828F10J		10FR	
GC05828F12J		12FR	
GC05828F14J		14FR	

NETELON CATHETER



Lateral Drilling

Product code	Description	Size	Tube Length
GC05802F10E-GC05802F16E	Female	FR10-FR16	23.5cm
GC05802M08B-GC05802M20B	Male	FR08-FR20	38.5cm
GC05802M06A-GC05802M20A	Male	FR06-FR20	42cm

Distal Drilling

Product code	Description	Size	Tube Length
GC05801F10C-GC05801F16C	Female	FR10-FR16	29cm
GC05801M06G-GC05801M14G	Male	FR06-FR14	40cm
GC05801M10D-GC05801M16D	Male	FR10-FR16	44cm

HYDROPHILIC CATHETER

Hydrophilic catheters are a type of pre-lubricated single-use, urinary catheter with a special coating applied to the catheter during the manufacturing process. When introduced to sterile water or sterile saline, the hydrophilic coating is activated and the catheter is ready for insertion.



Product code	Description
GC05805M	Male, 08-18FR, with water sachet
GC05805F	Female, 08-18FR, with water sachet

CATHETER SPIGOT

Catheter spigot is used to provide stop of flow for catheters hygienically during nursing procedure. It is non-invasive utilized to seal the catheter for a short time to allow urine collect in bladder.

FEATURES:

- Graduated design makes it easily accommodate tubing of various internal diameters
- Ribbed design is available
- Two sizes are available for different tube diameter




Product code	Spec
GC101301	Ribbed type
GC1071	O.D. 10mm
GC1072	O.D. 13mm
GC1939A01	Foley catheter adaptor

RECTAL CATHETER


The rectal plastic catheters are used for examination of rectum, administration of medicines and rectal drainage. One end of the catheter has central opening and the other has color plug. Rectal catheters are used in conjunction with delivery systems to provide an effective method of keeping your bowel movements regular and predictable.

- 04FR Grey
- 06FR Light green
- 08FR Light blue
- 10FR Black
- 12FR White
- 14FR Green
- 16FR Orange
- 18FR Red
- 20FR Yellow
- 22FR Purple
- 24FR Deep blue



Product code	Description	Size	Tip Type
GC050A04-GC050A24	Connector with cap	04-24FR	Open tip, with one smooth distal lateral eye
GC050B26-GC050B32	Without connector	26-32FR	

Product code	Size	Description
GC050C06	06FR	Funnel connector, open tip, one lateral eye, 06FR
GC050C08	08FR	Funnel connector, open tip, one lateral eye, 08FR
GC050C10	10FR	Funnel connector, open tip, one lateral eye, 10FR
GC050C12	12FR	Funnel connector, open tip, one lateral eye, 12FR
GC050C14	14FR	Funnel connector, open tip, one lateral eye, 14FR
GC050C16	16FR	Funnel connector, open tip, one lateral eye, 16FR
GC050C18	18FR	Funnel connector, open tip, one lateral eye, 18FR
GC050C20	20FR	Funnel connector, open tip, one lateral eye, 20FR



- 18FR Red
- 20FR Yellow
- 22FR Purple
- 24FR Deep blue
- 26FR White
- 28FR Grey
- 30FR Black

Product code	Description	Size	Tip Type
GC050D18-GC050D30	Funnel connector	18-30FR	Open tip, with one smooth distal lateral eye



SILICONE CATHETER

The Silicone catheter is a two-way catheter indicated for routine drainage of the bladder. It is a urinary catheter intended for short as well as long term drainage of the bladder. For routine and postoperative drainage.



Product code	Description
GC062A0601	06FR, 3ml, length 30cm, with guide-wire
GC062A0621	06FR, 3ml, length 55cm with guide-wire
GC062A0801	08FR, 3ml, length 30cm with guide-wire
GC062A0811	08FR, 3ml, length 35cm with guide-wire
GC062A0822	08FR, 5ml, length 55cm with guide-wire
GC062A1001	10FR, 3ml, length 30cm with guide-wire
GC062A1011	10FR, 3ml, length 35cm with guide-wire
GC062A1022	10FR, 5ml, length 55cm with guide-wire
GC062A1213	12FR, 10ml, length 35cm
GC062A1223	12FR, 10ml, length 55cm

SILICONE MALE CATHETER

Silicone male catheter is a self-adhering catheter that is made from odor-free silicone. It is latex-free to reduce skin irritation.

Transparent and clear, allowing visual inspections without the need to remove the catheter. Built in adhesive makes applications easier and allows longer wear time. A kink-proof funnel is designed to facilitate uninterrupted urine flow.



Size Indicator



Hair Protector

Product code	Description
GC063A25	I.D. 25mm
GC063A29	I.D. 29mm
GC063A32	I.D. 32mm
GC063A36	I.D. 36mm
GC063A41	I.D. 41mm
GC063B25	I.D. 25mm. With hair protector and size indicator
GC063B29	I.D. 29mm. With hair protector and size indicator
GC063B32	I.D. 32mm. With hair protector and size indicator
GC063B36	I.D. 36mm. With hair protector and size indicator
GC063B41	I.D. 41mm. With hair protector and size indicator

03




IRRIGATION SET

√ BLADDER IRRIGATION SET16-21
MAPLE Irrigation Set16
PEONY Irrigation Set16
SANDALWOOD Irrigation Set17
FREESIA Irrigation Set17
AZALEA Irrigation Set18
BAMBOO Irrigation Set18
JASMINE Irrigation Set19
HOLLY Irrigation Set19
DAFFODIL Irrigation Set20-21



BLADDER IRRIGATION SET

MAPLE Irrigation Set



Product code	GC081101	GC081102	GC081104
Non vented spike	1	2	4
Drip chamber	1	1	1
Click clamp		2	4
Roller regulator	1	1	1
Soft tube adaptor	1	1	1
Drawing			





PEONY Irrigation Set

Product code	GC081201	GC081202
Non vented spike	1	2
Pillar drip chamber	1	1
Click clamp		2
Roller regulator	1	1
Soft tube adaptor	1	1
Drawing		




SANDALWOOD Irrigation Set



Product code	GC081501	GC081511	GC081502	GC081512
Non vented spike	1	1	2	2
Pillar drip chamber	1	1	1	1
Click clamp			2	2
Roller regulator	1	1	1	1
Soft tube adaptor	1	1	1	1
Luer lock connector		1		1
Drawing				

FRESIA Irrigation Set

Product code	GC081521
Non vented spike	1
Drip chamber	1
Roller regulator	1
Soft tube adaptor	1
Drawing	



AZALEA Irrigation Set



Product code	GC081401	GC081402	GC081404
Non vented spike	1	2	4
Pillar drip chamber	1	1	1
Click clamp		2	4
Roller regulator	1	1	1
Soft tube adaptor	1	1	1

Drawing			
---------	--	--	--

BAMBOO Irrigation Set



Product code	GC082302	GC082304
Non vented spike	2	4
Jumbo drip chamber	1	1
Click clamp	2	4
Roller regulator	1	1
Soft tube adaptor	1	1

Drawing		
---------	--	--





JASMINE Irrigation Set



Product code	GC081702	GC081722
Non vented spike	2	2
Drip chamber	1	1
Roller regulator	1	1
Soft tube adaptor	1	1
Luer lock connector		1
Drawing		





HOLLY Irrigation Set

Product code	GC081801	GC081802
Non vented spike	1	2
Drip chamber	1	1
Roller regulator	1	1
Soft tube adaptor	1	1
Drawing		



DAFFODIL Irrigation Set

Product code	Description	Picture
GC081601	One-way, drip chamber with safety lock, 1.85M tubing with roller regulator and connector	
GC081601	One-way, drip chamber with safety lock, 3.2M tubing with roller regulator and connector	
GC081612	Two-way, with roller regulators, 1.2M tubing with drip chamber and connector	

Product code	Description	Picture
GC081621	One-way, drip chamber with safety lock, 3.2M tubing with roller regulator, connected with Luer lock connector with Luer cap	
GC081631	One-way, non-vented spike with drip chamber, 3.2M tubing with roller regulator, connected with Luer lock connector and Luer cap	
GC081641	Regulator and pillar drip chamber, connected with connector	
GC1964E01	Endoscopy irrigation set with Luer lock connector and Luer cap	

04

**SINGLE-USE DIGITAL
FLEXIBLE URETEROSCOPE**

√ SINGLE-USE DIGITAL FLEXIBLE URETEROSCOPE22

SINGLE-USE DIGITAL FLEXIBLE URETEROSCOPE

FEATURES:

Endoscopic examination or therapeutic surgery for the ureter and renal pelvis. It delivers detailed high-resolution digital images for high-quality visualization and navigation - to help you remove stones quickly, safely and affordably.



Product code	Working length	Bending angle	Pixel camera	Maximum outer diameter	Working channel	Coated or not
GC090B01	650mm	U275 °/D275 °	160000	3.1mm	1.2mm	
GC090B02	650mm	U275 °/D275 °	160000	2.8mm	1.2mm	
GC090B03	650mm	U275 °/D275 °	160000	2.5mm	1.2mm	
GC090H01	650mm	U275 °/D275 °	160000	3.1mm	1.2mm	✓
GC090H02	650mm	U275 °/D275 °	160000	2.8mm	1.2mm	✓
GC090H03	650mm	U275 °/D275 °	160000	2.5mm	1.2mm	✓



CONTACT US



**LOUJIN INDUSTRIAL PARK, SHUOFANG,
WUXI 214143, Jiangsu, P.R.CHINA**



86-510-82733575



INFO@GCMEDICA.COM



WWW.GCMEDICA.COM

Version: DJA23

We reserve the right to correct any errors that may occur in this brochure.