




Adjustable Tracheostomy Tube

- ✔ The fast locking mechanism on flange easily secures tube to proper depth.
- ✔ Supplied with Tapered obturator aids quick insertion.

Adjustable Tracheostomy Tube



Adjustable flange can facilitate better fit to the patient's anatomy.

- ◆ The fast locking mechanism on flange easily secures tube to proper depth.
- ◆ Supplied with Tapered obturator aids quick insertion.
- Soft & flexible flange provided with inner & outer diameter sizes
- Standard 15mm connector conforms to ISO standards, can rotate 360 degree & enable access to the anesthetic circuit from any direction
- Supplied with cotton neck tape

No.	Product Code - Uncuffed	Product Code - Cuffed	Size
01	GC8041C030	GC8041D030	3.0
02	GC8041C035	GC8041D035	3.5
03	GC8041C040	GC8041D040	4.0
04	GC8041C045	GC8041D045	4.5
05	GC8041C050	GC8041D050	5.0
06	GC8041C055	GC8041D055	5.5
07	GC8041C060	GC8041D060	6.0
08	GC8041C065	GC8041D065	6.5
09	GC8041C070	GC8041D070	7.0
10	GC8041C075	GC8041D075	7.5
11	GC8041C080	GC8041D080	8.0
12	GC8041C085	GC8041D085	8.5
13	GC8041C090	GC8041D090	9.0
14	GC8041A095	GC8041D095	9.5
15	GC8041A100	GC8041D100	10.0