



ENT

SINGLE USE SUCTION TUBE

Product
introduction



CONTENTS

01	Polished Suctions	01
02	Polished Suctions with Obturator	02
03	Rounded Suctions	04
04	Fenestrated Suctions	05
05	Rounded & Fenestrated Suctions	06
06	Rounded Suction with 4 Fenestrations	07
07	Soft Tip Suction	08
08	Dissector Suction	08



ENT SINGLE USE SUCTION TUBE

ENT Single Use Suction Tube are used in a variety of medical procedures including laryngeal suction, where they irrigate and aspirate the throat, and aural micro-suction, where they are used to remove stubborn or impacted earwax. A range of supporting accessories allow for greater clinical control.


Our core range of ENT Single Use Suction Tube are available in a variety of tube lengths and lumen sizes to suit each case. Selected suction tubes are available with atraumatic tips, with a choice of polished, rounded, fenestrated, rounded tips and so on. Our suction tubes also available for Zoellner, Mastoid, Frazier, Magill.

The Zoellner suction tubes are used in aural micro-suction to remove impacted earwax.

The Mastoid suction tubes are used in radical procedures on the frontal sinuses, in operations on the antrum, in external ethmoidectomy and also in evisceration of the eye and operations on the lacrimal.

The Frazier suction tubes are used to remove fluid or debris from nasal cavities and in laryngeal suction.

The Magill suction tubes are used to remove impacted earwax during aural micro-suction.



Polished Tip

Polished Suctions

All suction tubes have smooth polished tips. Ergonomic balanced design for ease of use. Widely range of lengths and lumens to suit each clinic use.

Purple Handle Is Available



F Features:

- Single use for ENT surgery.
- A variety of tube lengths and lumen sizes for each clinic suit.
- Lightweight, non-blocking design.
- Ergonomic design polymer handle comfortably suits left or right-handed users.
- Available with control hole handle, allows greater suction when covered with user's finger.

Product Code	Size	Description
1262S20W1	Working Length 85/60mm, OD=2mm	With Control vent, Zoellner
1262S30W1	Working Length 85/60mm, OD=3mm	With Control vent, Mastoid
1262S20W2	Working Length 128/90mm, OD=2mm	With Control vent, Frazier
1262S30W2	Working Length 128/90mm, OD=3mm	With Control vent, Frazier


1262S40W2	Working Length 128/90mm, OD=4mm	With Control vent, Frazier
1262S50W2	Working Length 128/90mm, OD=5mm	With Control vent, Frazier
1262S40W3	Working Length 138/100mm, OD=4mm	With Control vent, Frazier
1262S20W4	Working Length 168/130mm, OD=2mm	With Control vent, Frazier
1262S30W4	Working Length 168/130mm, OD=3mm	With Control vent, Frazier
1262S40W4	Working Length 168/130mm, OD=4mm	With Control vent, Frazier
1262S50W4	Working Length 168/130mm, OD=5mm	With Control vent, Frazier
1262R20W0	Working Length 85/35mm, OD=2mm	With Control vent, Magill
1262R30W0	Working Length 85/35mm, OD=3mm	With Control vent, Magill
1264R20W0	Working Length 85/35mm, OD=2mm	Purple handle With Polished Tip, Frazier
1264S20W1	Working Length 85/60mm, OD=2mm	Purple handle With Polished Tip, Frazier
1264S30W2	Working Length 128/90mm, OD=3mm	Purple handle With Polished Tip, Frazier



Polished Suctions with Obturator

F Features:

- Single use for ENT surgery.
- Flexible cannula could be bent to needed position without occlusion.
- Available with integrated transparent handle for visual check.



Polished Tip with obturator

Product Code	Size	Description
1242A08	Fr08	With flexible cannula and obturator, Frazier
1242A10	Fr10	With flexible cannula and obturator, Frazier
1242A12	Fr12	With flexible cannula and obturator, Frazier
1242A18	Fr18	With flexible cannula and obturator, Frazier

Transparent Handle Polished Suctions With Obturator



Product Code	Size	Description
1241A06	Fr06	Transparent handle, With flexible cannula and obturator, Frazier
1241A08	Fr08	Transparent handle, With flexible cannula and obturator, Frazier
1241A10	Fr10	Transparent handle, With flexible cannula and obturator, Frazier
1241A12	Fr12	Transparent handle, With flexible cannula and obturator, Frazier
1241B06	Fr06	Transparent handle, With flexible cannula and X-ray line obturator, Frazier
1241B08	Fr08	Transparent handle, With flexible cannula and X-ray line obturator, Frazier
1241B10	Fr10	Transparent handle, With flexible cannula and X-ray line obturator, Frazier
1241B12	Fr12	Transparent handle, With flexible cannula and X-ray line obturator, Frazier



Rounded Suctions

A range of suction tubes with a smooth rounded tip for atraumatic tissue contact. Used for ENT (Ear-Nose-Throat) surgery.



F Features:

- Fine smooth rounded tip for atraumatic tissue contact.
- Available with a variety of lengths and lumens for each cases.




Rounded Tip

Product Code	Size	Description
1262P40W1	Working Length 85/60mm, OD=4mm	With Rounded Tip, Frazier
1262P20W2	Working Length 128/90mm, OD=2mm	With Rounded Tip, Frazier
1262P30W2	Working Length 128/90mm, OD=3mm	With Rounded Tip, Frazier
1262P40W2	Working Length 128/90mm, OD=4mm	With Rounded Tip, Frazier
1262P20W4	Working Length 168/130mm, OD=2mm	With Rounded Tip, Frazier
1262P30W4	Working Length 168/130mm, OD=3mm	With Rounded Tip, Frazier
1262P40W4	Working Length 168/130mm, OD=4mm	With Rounded Tip, Frazier
1262O30W0	Working Length 85/35mm, OD=3mm	With Rounded Tip, Magill



Fenestrated Suctions

A range of suction with fenestrated tips to disperse aspiration and prevent suction trauma.



Fenestrated Tip

Product Code	Size	Description
1262T20W1	Working Length 85/60mm, OD=2mm	With Fenestrated Tip, Frazier
1262T30W1	Working Length 85/60mm, OD=3mm	With Fenestrated Tip, Frazier
1262T40W1	Working Length 85/60mm, OD=4mm	With Fenestrated Tip, Frazier
1262T20W2	Working Length 128/90mm, OD=2mm	With Fenestrated Tip, Frazier
1262T30W2	Working Length 128/90mm, OD=3mm	With Fenestrated Tip, Frazier
1262T20W4	Working Length 168/130mm, OD=2mm	With Fenestrated Tip, Frazier
1262T30W4	Working Length 168/130mm, OD=3mm	With Fenestrated Tip, Frazier
1262T40W4	Working Length 168/130mm, OD=4mm	With Fenestrated Tip, Frazier



Rounded & Fenestrated Suctions

Suctions that benefit from a tip that is both fenestrated and rounded.

F Features:

- Suction control.
- For micro-laryngoscopy (MLS).
- Rounded opening for suction protection.



Rounded & Fenestrated Tip

Product Code	Size	Description
1262Q20W2	Working Length 128/90mm, OD=2mm	With Rounded & Fenestrated Tip, Frazier
1262Q30W2	Working Length 128/90mm, OD=3mm	With Rounded & Fenestrated Tip, Frazier
1262Q40W2	Working Length 128/90mm, OD=4mm	With Rounded & Fenestrated Tip, Frazier
1262Q20W4	Working Length 168/130mm, OD=2mm	With Rounded & Fenestrated Tip, Frazier
1262Q30W4	Working Length 168/130mm, OD=3mm	With Rounded & Fenestrated Tip, Frazier
1262Q40W4	Working Length 168/130mm, OD=4mm	With Rounded & Fenestrated Tip, Frazier




Fenestrated Suctions

This ultra-delicate suction has 4 fenestrations that are 1mm in diameter at the tip.



Rounded Suction with 4 Fenestrations Tip

Product Code	Size	Description
1262H40W2	Working Length 128/90mm, OD=4mm	Rounded Tip, Fenestrations 4*1 mm, Mastoid



Soft Tip



Soft Tip Suction

A suction with a soft silicone tip for very delicate procedures.

Product Code	Size	Description
1262M20W2	Working Length 128/90mm, OD=2mm	With Soft Tip



Dissector Suction

A straight tube finished with a flattened polished tip.



Dissector Flattened Tip



Product Code	Size	Description
1262N40W4	Working Length 168/130mm, OD=4mm	With Dissector Flattened Tip



ENT

SINGLE USE SUCTION
TUBE

CONTACT US



0086-510-88880366



LOUJIN INDUSTRIAL PARK, SHUOFANG,
WUXI, 214143, P.R.CHINA



INFO@GCMEDICA.COM



WWW.GCMEDICA.COM
