


UROLOGY & GYNAECOLOGY



Adress: Loujin Industrial park, Shuofang, Wuxi, 214143, P.R.China



Tel: 86 510 8273 3513 Fax: 86 510 8273 0717



E-mail: info@gcmedica.com



Website: www.gcmedica.com



GCMEDICA ENTERPRISE LTD.,(WUXI)



GCMEDICA PROFILE

GCMEDICA Group founded in 1998 with a mission dedicating to all medical personnel and patients, combines many years of experience and expertise, GCMEDICA specializes in the manufacturing scientific research and development of disposable medical devices and advanced moulds. GCMEDICA have built up a leading group company concentrated in various products production with CE certificate and registered with FDA US. With persistent and diligent exploration, GCMEDICA extents the product ranges to Respiratory, Anesthesia, Urology, Surgery, Operating room, Gynecology products.

GCMEDICA located in new industrial zone-Wuxi, which owns excellent convenient transportation. GCMEDICA has building area of 10000m² manufacture base which incorporates a 100,000-class clean room of 5500m² and well-equipped laboratories. GCMEDICA mainly manufactures Class I, Class II and Class III disposable medical devices. The applied expertise of our R&D team ensures our latest products are always innovative.

GCMEDICA is a professional and 100% purely manufacturer, who with the whole manufacturing process completed within our plants, from design, manufacture to sterilization.

GCMEDICA upholds “devoted to improve the lives of the people” as its spirit of enterprise. Our commitment to social engagement and responsibility is an important aspect of our company culture and philosophy.

GCMEDICA keeps the concept of being oriented by science, technology and taking quality as life. Providing the reliable products and credit has always developed its advantage. Making the products as a faith, observing market is guides to lead GCMEDICA a world famous brand.




CONTENTS

◆ UROLOGICAL	1~13
※ CYSTO/BLADDER IRRIGATION SET	1~5
• M-Easy Bladder Irrigation Set	1
• B-Cylind Bladder Irrigation Set	1
• S-Tur Bladder Irrigation Set	2
• S-Uni Bladder Irrigation Set	2
• B-Uro Bladder Irrigation Set	3
• Premi Bladder Irrigation Set	3
• J-Pump Bladder Irrigation Set	4
• J-Tur Bladder Irrigation Set	4
• H-Pump Bladder Irrigation Set	5
• Sup-Flow Bladder Irrigation Set	5
※ ARTHROSCOPY IRRIGATION SET	6
※ DRAINAGE BAG WITH URINE METER	7
※ LOW BED URINE BAG	7
※ URINARY DRAINAGE BAG.....	8
※ ENEMA KITS	9
※ SITZ BATH KITS	9
※ URINARY DRAINAGE LEG BAG	10
※ CLICK SEAL SPECIMEN CONTAINER	11
※ SUCTION SPECIMEN DRAIN BAG	11
※ NELATON CATHETER	12
※ SPIGOT CATHETER AND ADAPTOR	12
※ SILICON MALE CATHETER	13
※ RECTAL BALLOON CATHETER	13
◆ GYNAECOLOGY	14~16
※ UTERINE ASPIRATION TUBING	14
※ UTERINE ASPIRATION CURETTE	14
※ VAGINAL SPECULUM	15~16
※ AMNIO HOOK	16
※ IUD REMOVAL HOOK	16

CYSTO / BLADDER IRRIGATION SET



M-Easy Bladder Irrigation Set



Product code	081101	081102	081104
Non vented spike	1	2	4
Pillar drip chamber	1	1	1
Click clamp		2	4
Y-connector		1	3
Roller regulator	1	1	1
Silicon tube adapter	1	1	1
Luer lock connector			
Drawing			





B-Cylind Bladder Irrigation Set




Product code	081201	081202
Non vented spike	1	2
Pillar drip chamber	1	1
Click clamp	2	
Y-connector		
Roller regulator	1	1
Silicon tube adapter		1
Drawing		

S-Tur Bladder Irrigation Set



Product code	081501	081511	081502	081512
Non vented spike	1	1	2	2
Drip chamber	1	1	1	1
Click clamp			2	2
Y-connector		1		1
Roller regulator	1	1	1	1
Silicon tube adapter	1	1	1	1
Lure lock connector		1		1
Drawing				

S-Uni Bladder Irrigation Set

Product code	081521
Non vented spike	1
Drip chamber	1
Roller regulator	1
Silicon tube adapter	1
Drawing	




B-Uro Bladder Irrigation Set



Product code	081401	081402	081404
Non vented spike	1	2	4
Drip chamber	1	1	1
Click clamp		2	4
Y-connector		1	3
Roller regulator	1	1	1
Silicon tube adapter	1	1	1

Drawing



Premi Bladder Irrigation Set

Product code	082302	082304
Non vented spike	2	4
Jumbo drip chamber	1	1
Click clamp		
Y-connector		3
Roller regulator		1
Silicon tube adapter	1	1

Drawing




J-Pump Bladder Irrigation Set



Product code	081702	081722
Non vented spike	2	2
Drip chamber	1	2
Roller regulator	1	1
Silicon tube adapter	1	1
Luer lock connector		1

Drawing



J-Tur Bladder Irrigation Set





Product code	081712	081732
Non vented spike	2	2
Roller regulator	1	1
Silicon tube adapter	1	1
Luer lock		1

Drawing



H-Pump Bladder Irrigation Set




Product code	081801	081802
Non vented spike	1	2
Drip chamber	1	1
Roller regulator	1	1
Silicon tube adapter	1	1
Luer lock connector		1
Drawing		

Sup-Flow Bladder Irrigation Set

FEATURE:





·Burette drip chamber and extra wider tube I.D. for large flow irrigation.



Product code	084404
Non vented spike	4
Drip chamber	1
Roller regulator	1
Silicon tube adapter	1
Drawing	

ARTHROSCOPY IRRIGATION SET



Product code	080601	080602	080611	080604
Non vented spike	2	2	2	4
Drip chamber	1	2	1	1
Click clamp	3	3	3	5
Y-connector	1	1	1	3
Roller clamp			1	
Silicon tube adapter	1	1	1	1
Drawing				

DRAINAGE BAG WITH URINE METER

Drainage Bag with Urine Meter is a closed system used to accurately measure urine output and collect the urine. The device features a rigid, clear and transparent measured volume chamber, and a 2600mL soft drainage bag.

Product code 051E01

FEATURES:

- 2600ml drainage bag with 500ml meter
- Anti-reflux valve
- Access port
- Double wrapped
- With convenient hanger
- 150cm tubing
- Latex free
- Single patient use



LOW BED URINE BAG

Low-Bed Drainage Bag features to facilitate proper drainage and anti-reflux valve for back-flow prevention. The exclusive slide-tap drainage port allows for minimal splashing during drainage. The drainage Bag has a reversible design for versatile placement on either side of the bed. This drain bag hangs only 6-inches (18cm) to allow low placement on an adjustable bed without touching the floor.



Product code 051M01

FEATURES:

- With cross valve. With anti-reflux valve
- With reinforced double hanger & rope hanger

URINARY DRAINAGE BAG

- 2000ml,with T-Tap outlet valve,and soft tube pocket
- Needless sample port with cap
- Two parts anti-reflux drip chamber,air-vent on the bag
- With reinforced double hanger & rope hanger
- With bed sheet clamp
- The tubing is \varnothing 10.0×1100mm
- Latex-free,packed sterile

Product code 051B01

FEATURES:

- Diamond shaped bag
- With cross valve. With anti-reflux valve
- With reinforced single hanger



Product code 051B02

FEATURES:

- Pear-shaped bag
- With reinforced single hanger & rope hanger

Product code 051D01

FEATURES:

- Triangle shaped bag
- With single reinforced hanger



ENEMA KITS

FEATURES:

- 1500ml bag or bucket with 60 inches tubing, with pre-lubricated tip
- Soap packet
- Water-proof drape: 420x420mm



Product code	Description
056A15	Enema bag only
056B15	Enema kit with bucket

SITZ BATH KITS

FEATURES:

- Ideal for the treatment of hemorrhoids or other perineal conditions.
- Provides wide, contoured edges for maximum comfort and support.
- Fits most standard toilets and commodes.
- Water bag can be positioned on toilet seat or on top of toilet tank.
- Easy to clean durable plastic basin.
- 2000ml bag with 60 inches tubing.
- Fits standard toilets and commodes.
- Portable and can administrate soothing fluids to perinea area.



Product code	Description
054A01	Sitz bag only
059A01	Sitz bag with a bath basin

URINARY DRAINAGE LEG BAG

Product code	Spec
055C06	600ml
055C16	600ml, with extension tube
055C09	900ml
055C19	900ml,with extension tube

FEATURES:

- Square shape,with anti-reflux valve
- Flip style valve drainage outlet
- Two velcro straps
- With twist valve
- With extension tubing



Product code	Spec
055D04	450ml
055D07	750ml
055D09	900ml

FEATURES:

- Square shape,with anti-reflux valve
- Twist style valve drainage outlet
- With double high elastic belts

Product code	Spec
055B05	500ml
055B07	700ml
055B10	1000ml

Urinary Leg Bag is durably constructed for odor containment and leak resistance. The anti-reflux chamber reduces gross reflux.

FEATURES:

- Elliptic shape , with anti-reflux valve
- Equipped with handle/hanger
- Design allows for even filling and draining
- Includes two flexible straps
- Single use
- Latex free



CLICK SEAL SPECIMEN CONTAINER

FEATURES:

- With a sound of click when cap is sealed tight, decreasing the risk of contamination.
- Special cap design ensures convenience and easy job for nurses.
- With graduation for fast volume reading.
- Available with 90ml and 120ml capacity.



Product code	Description
061A03	90ml
061A04	120ml

SUCTION SPECIMEN DRAIN BAG


FEATURES:

- Open mesh filter allow for fast drainage of fluids
- Durable, latex drape has flexible hose
- Velcro securement eliminates the need to use tape



Product code	05701	05702
Filter	one-piece filter	removable filter
Flexible hose	with	with

Drawing



NELATON CATHETER



Product code	Description
05801F	26cm length tube with smooth two eyes, funnel connector
05802F	26cm length tube with two eyes, funnel connector
05801M	40cm length tube with smooth two eyes, funnel connector
05802M	40cm length tube with two eyes, funnel connector
058A	40cm length tube with one eye, closed tip, connector with cap

SPIGOT CATHETER AND ADAPTOR

Catheter Spigot is used to provide stop of flow for catheters hygienically during nursing procedure. It is non-invasive utilized to seal the catheter for a short time to allow urine collect in bladder.

FEATURES:

- Graduated design makes it easily accommodate tubing of various internal diameters.
- Ribbed design is available.
- Two sizes are available for different tube diameter.




Product code	Description
1971	0.D.10mm
1972	0.D.13mm
101301	Ribbed type
1939A01	Foley catheter adaptor

SILICON MALE CATHETER

FEATURES:

Silicon male catheter is a self-adhering catheter that is made from odor-free silicone. It is latex-free to reduce skin irritation. Transparent and clear, allowing visual inspections without the need to remove the catheter. Built in adhesive makes applications easier and allows longer wear time. A kink-proof funnel is designed to facilitate uninterrupted urine flow.

- Self-Adhering Bond.
- Kink Proof Funnel for Uninterrupted Flow.
- Latex-Free to Avoid Irritation.
- Constructed with 100% Silicone.
- Transparent for Easy Visual Skin Inspections.
- Hair protector and size indicator is available for convenient application.




Product code	Description
063A25	I.D. 25mm
063A29	I.D. 29mm
063A32	I.D. 32mm
063A36	I.D. 36mm
063A41	I.D. 41mm
063B25	I.D. 25mm, With hair protector and size indicator
063B29	I.D. 29mm, With hair protector and size indicator
063B32	I.D. 32mm, With hair protector and size indicator
063B36	I.D. 36mm, With hair protector and size indicator
063B41	I.D. 41mm, With hair protector and size indicator

RECTAL BALLOON CATHETER

FEATURES:

- Used for measuring rectal pressure during urodynamic studies.
- Silicone construction avoids latex sensitivity concerns.
- Made of PVC with silicone balloon, single lumen
- Universal adapter with sealing cap
- Sterile, single use




Product code	6035C01
--------------	---------

UTERINE ASPIRATION TUBING

FEATURES:

- Flexible cuff for convenient operations.
- Special wide bore tubing, thick walled, high vacuum PVC tubing-proof seals.
- Available with flexible cuff, providing extra flexibility at crucial point.
- Easily connects to Curettes.




Product code	Description
7050A01	Rotary handle only
7051A01	I.D.12mm bore with rotary handle
7051B01	I.D.12mm with rotary handle and adapter
7051C01	I.D.12mm with rotary handle, adapter and flexible cuff

UTERINE ASPIRATION CURETTE

FEATURES:

- Rigid, clear, plastic curettes for uterine aspiration
- Molded markings for size identification
- Available with straight one and curved one
- Available sizes: 6mm, 7mm, 8mm, 9mm, 10mm, 11mm, 12mm, 14mm, 16mm



Product code	Description
7049A	Straight
7049B	Curved

VAGINAL SPECULUM

FEATURES:


- Crystal clear blade for accurate gynecologic checking
- Unbreakable design for patient's safety.
- Smooth and rounded end for patient safety and comfort.
- Transparent finish yields excellent visibility, sampling and examination accuracy.
- With middle screw type design.



Product code	Description
7048L01	Size L
7048L02	Size L, unbreakable
7048M01	Size M
7048M02	Size M, unbreakable
7048S01	Size S
7048S02	Size S, unbreakable

FEATURES:

- Made from non-toxic material.
- Available with common material and unbreakable.
- Middle screw type.



Product code	Description
7047L01	Size L, french type
7047M01	Size M, french type
7047S01	Size S, french type
7047XS01	Size XS, french type

FEATURES:

- Made from non-toxic material.
- Available with common material and unbreakable.
- Middle screw type.



Product code	Description
7045L01	Size L, side screw type
7045M01	Size M, side screw type
7045S01	Size S, side screw type

AMNIO HOOK

FEATURES:

- Used during vaginal deliveries to puncture the amniotic sac.
- Disposable to prevent infection.



Product code	Description
7053A01	

IUD REMOVAL HOOK

FEATURES:

- For easy and effective retrieval and removal of an IUD.
- Disposable to prevent infection.



Product code	Description
7057A01	